

CAUTIONS

U.S. Federal Law restricts sale and use to or on order of a Physician.

Gamma Sterilization is performed on the CleanTouch pouches. Never open the CleanTouch pouch until just before use. Do not use the CleanTouch product if the outer pouch has been broken. The sterilization effect may be weakened and the product may become non-sterile.

WARNINGS

This product is designed for single use. Do not re-use. Re-use will result in sub-optimal performance of the PVA material and may contaminate the surgical field.

If the expiration date has passed, do not use. Use past the expiration date may result in sub-optimal performance of the PVA material and the product may become non-sterile.

© 2020 Copyright MED-LOGICS®, Inc.
All rights reserved.

P/N: 900-0795-011

CLEANTOUCH
PVA EYE SPONGES



ML100

Clean-Compressed PVA Eye Spears

 **MED-LOGICS, INC.**

CE
2797

R
ONLY

CLEANTOUCH
PVA EYE SPONGES



ML100

Clean-Compressed PVA Eye Spears

 **MED-LOGICS, INC.**

CE
2797

R
ONLY

CAUTIONS

U.S. Federal Law restricts sale and use to or on order of a Physician.

Gamma Sterilization is performed on the CleanTouch pouches. Never open the CleanTouch pouch until just before use. Do not use the CleanTouch product if the outer pouch has been broken. The sterilization effect may be weakened and the product may become non-sterile.

WARNINGS

This product is designed for single use. Do not re-use. Re-use will result in sub-optimal performance of the PVA material and may contaminate the surgical field.

If the expiration date has passed, do not use. Use past the expiration date may result in sub-optimal performance of the PVA material and the product may become non-sterile.

© 2020 Copyright MED-LOGICS®, Inc.
All rights reserved.

P/N: 900-0795-011

Description of Product

The CleanTouch PVA Eye Spears are sterile, single-use devices designed to absorb excess fluids on the human eye.



MED-LOGICS, Inc.®
1627 Enterprise Street, Athens, TX 75751 USA
Tel: +1-949-582-3891/ E-mail: info@mlogics.com



European Representative
E C Rep Ltd,
5 Fitzwilliam Square East, Dublin 2, D02 R744, Ireland
Tel: +353-89 225-1951/ E-mail: info@ecrep.ie



UK Representative
Responsible Person Ltd,
Summit House, 4-5 Mitchell Street, Edinburgh, EH6 7BD, UK
Tel: +44-7543-672-888/ E-mail: info@responsible-person.co.uk

Instructions for Use

1. Remove foil pouches from the sterile pouch.
2. Remove spears from the foil pouch.
3. Introduce PVA material to any excess fluids. The PVA material will absorb fluid until it reaches its capacity.
4. Once a spear reaches its capacity, discard and use a new spear.

Storage Conditions

The CleanTouch PVA Eye Spears have no special storage conditions, but should be stored under normal warehouse conditions, which includes protection from moisture, extreme cold (less than 5°C), and excessive heat (greater than 40°C).

Description of Product

The CleanTouch PVA Eye Spears are sterile, single-use devices designed to absorb excess fluids on the human eye.



MED-LOGICS, Inc.®
1627 Enterprise Street, Athens, TX 75751 USA
Tel: +1-949-582-3891/ E-mail: info@mlogics.com



European Representative
E C Rep Ltd,
5 Fitzwilliam Square East, Dublin 2, D02 R744, Ireland
Tel: +353-89 225-1951/ E-mail: info@ecrep.ie



UK Representative
Responsible Person Ltd,
Summit House, 4-5 Mitchell Street, Edinburgh, EH6 7BD, UK
Tel: +44-7543-672-888/ E-mail: info@responsible-person.co.uk

Instructions for Use

5. Remove foil pouches from the sterile pouch.
6. Remove spears from the foil pouch.
7. Introduce PVA material to any excess fluids. The PVA material will absorb fluid until it reaches its capacity.
8. Once a spear reaches its capacity, discard and use a new spear.

Storage Conditions

The CleanTouch PVA Eye Spears have no special storage conditions, but should be stored under normal warehouse conditions, which includes protection from moisture, extreme cold (less than 5°C), and excessive heat (greater than 40°C).

Symbols Used on Labeling



Catalog Number



Lot Number



Manufacturer



Authorized Representative in the EU



Authorized Representative in the UK



CE Marking and Authorized Rep No.



Consult Instructions for Use



Caution



Do Not Reuse



Do not use if package is damaged



Sterilized using Ethylene Oxide



For Use by, or on the Order of, a Physician



Date of Manufacture: (YYYY-MM)



Use by: (YYYY-MM)



Medical Device



Unique Device Identifier

Symbols Used on Labeling



Catalog Number



Lot Number



Manufacturer



Authorized Representative in the EU



Authorized Representative in the UK



CE Marking and Authorized Rep No.



Consult Instructions for Use



Caution



Do Not Reuse



Do not use if package is damaged



Sterilized using Ethylene Oxide



For Use by, or on the Order of, a Physician



Date of Manufacture: (YYYY-MM)



Use by: (YYYY-MM)



Medical Device



Unique Device Identifier